

Antibodies Product Datasheet

Product Name:	Phosphoserine Monoclonal Antibody(5B12)
Catalog Number:	N2027
Reactivity:	Species independent
Applications:	WB, IHC, IF, IP
Dilutions:	WB 1:500 - 1:2000, IHC 1:100 - 1:300, ELISA: 1:40000, IP 1:50-200, IF 1:50-200
Modification:	Unmodified/Total
Source:	Monoclonal, Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified Protein
Concentration:	
Gene Name:	
Protein Name:	Phosphoserine
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Phosphoserine protein.
Storage and Stability:	Store at -20° C/1 year