

Antibodies Product Datasheet

Product Name:	CaMKII $\beta / \gamma / \delta$ (Phospho Thr287) Monoclonal Antibody(4H2)
Catalog Number:	N2029
Reactivity:	Human, Rat, Mouse
Applications:	IHC, IF
Dilutions:	IHC 1:100-200. IF 1:50-200
Modification:	Unmodified/Total
Source:	Monoclonal, Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CaMKII $\beta / \gamma / \delta$ (Phospho Thr287)
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	CaMKII $\beta / \gamma / \delta$
Other Name:	CAMK2B; CAM2; CAMK2; CAMKB; Calcium/calmodulin-dependent protein kinase type II subunit beta; CaM kinase II subunit beta; CaMK-II subunit beta; CAMK2G; CAMK; CAMK-II; CAMKG; Calcium/calmodulin-dependent protein kinase type II subunit gamma;
Human Swiss Prot No.:	Q13554 / Q13555 / Q13557
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CaMKII $\beta / \gamma / \delta$ (Phospho Thr287) protein detects endogenous levels of CaMKII $\beta / \gamma / \delta$ (Phospho Thr287)
Storage and Stability:	Store at -20° C/1 year