

Antibodies Product Datasheet

Product Name:	Phospho GSK3 β (Ser9) Mouse mAb(6H6)
Catalog Number:	N2035
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	IHC1:50-200 ,WB 1:1000-2000
Modification:	Phospho Specific
Source:	Monoclonal, Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Phospho-GSK3 β (Ser9)
Concentration:	1 mg/ml
Gene Name:	GSK3B
Protein Name:	GSK3 β
Other Name:	Glycogen synthase kinase-3 beta (GSK-3 beta) (EC 2.7.11.26) (Serine/threonine-protein kinase GSK3B) (EC 2.7.11.1)
Human Swiss Prot No.:	P49841
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Phospho-GSK3 β (Ser9) at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year