

## Antibodies Product Datasheet

|                               |                                                                                                                                            |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | HA-Tag Polyclonal Antibody, HRP Conjugated                                                                                                 |
| <b>Catalog Number:</b>        | N2039                                                                                                                                      |
| <b>Reactivity:</b>            | Species independent                                                                                                                        |
| <b>Applications:</b>          | WB, ELISA                                                                                                                                  |
| <b>Dilutions:</b>             | Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:5000. |
| <b>Modification:</b>          | Unmodified                                                                                                                                 |
| <b>Source:</b>                | Polyclonal, Rabbit, IgG                                                                                                                    |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.                              |
| <b>Immunogen:</b>             |                                                                                                                                            |
| <b>Concentration:</b>         | 1mg/ml                                                                                                                                     |
| <b>Gene Name:</b>             | HA Tag                                                                                                                                     |
| <b>Protein Name:</b>          | HA-Tag                                                                                                                                     |
| <b>Other Name:</b>            | HA                                                                                                                                         |
| <b>Human Swiss Prot No.:</b>  |                                                                                                                                            |
| <b>Formulation:</b>           | Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.                                                     |
| <b>Specificity:</b>           | HA-Tag Polyclonal Antibody HRP Conjugated, specially designed for your Western blot analysis.                                              |
| <b>Storage and Stability:</b> | Store at -20° C/1 year                                                                                                                     |