

Antibodies Product Datasheet

Product Name:	Erbin (Phospho-Tyr1104) Antibody
Catalog Number:	P12-1021
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human Erbin (phospho Tyr1104)
Concentration:	1 mg/ml
Gene Name:	ERBB2IP
Protein Name:	Protein LAP2
Other Name:	ERBB2IP; ERBIN; KIAA1225; LAP2; Protein LAP2; Densin-180-like protein; Erbb2-interacting protein; Erbin
Human Swiss Prot No.:	Q96RT1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Erbin (Y1104) Polyclonal Antibody detects endogenous levels of Erbin protein only when phosphorylated at Y1104.
Storage and Stability:	-20°C/1 year