

Antibodies Product Datasheet

Product Name:	PDLIM5 (Phospho-Tyr251) Antibody
Catalog Number:	P12-1042
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human PDLIM5 (phospho Tyr251)
Concentration:	1 mg/ml
Gene Name:	PDLIM5
Protein Name:	PDZ and LIM domain protein 5
Other Name:	PDLIM5; ENH; L9; PDZ and LIM domain protein 5; Enigma homolog; Enigma-like PDZ and LIM domains protein
Human Swiss Prot No.:	Q96HC4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-PDLIM5 (Y251) Polyclonal Antibody detects endogenous levels of PDLIM5 protein only when phosphorylated at Y251.
Storage and Stability:	-20°C/1 year