

Antibodies Product Datasheet

Product Name:	Shb (Phospho-Tyr246) Antibody
Catalog Number:	P12-1055
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human Shb (phospho Tyr246)
Concentration:	1 mg/ml
Gene Name:	SHB
Protein Name:	SH2 domain-containing adapter protein B
Other Name:	SHB; SH2 domain-containing adapter protein B
Human Swiss Prot No.:	Q15464
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Shb (Y246) Polyclonal Antibody detects endogenous levels of Shb protein only when phosphorylated at Y246.
Storage and Stability:	-20°C/1 year