

Antibodies Product Datasheet

Product Name:	AQP3 Antibody
Catalog Number:	R12-2027
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AQP3. AA range:165-214
Concentration:	1 mg/ml
Gene Name:	AQP3
Protein Name:	Aquaporin-3
Other Name:	AQP3; Aquaporin-3; AQP-3; Aquaglyceroporin-3
Human Swiss Prot No.:	Q92482
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AQP3 Polyclonal Antibody detects endogenous levels of AQP3 protein.
Storage and Stability:	-20°C/1 year