

Antibodies Product Datasheet

Product Name:	ATP5J2 Antibody
Catalog Number:	R12-2032
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATP5J2. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	ATP5J2
Protein Name:	ATP5J2
Other Name:	ATP synthase f chain mitochondrial; ATP5JL; ATPK
Human Swiss Prot No.:	A6ND55
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATP5J2 Polyclonal Antibody detects endogenous levels of ATP5J2 protein.
Storage and Stability:	-20°C/1 year