

Antibodies Product Datasheet

Product Name:	BRCAA1 Antibody
Catalog Number:	R12-2043
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BRCAA1. AA range:761-810
Concentration:	1 mg/ml
Gene Name:	ARID4B
Protein Name:	AT-rich interactive domain-containing protein 4B
Other Name:	ARID4B; BRCAA1; RBBP1L1; RBP1L1; SAP180; AT-rich interactive domain-containing protein 4B; ARID domain-containing protein 4B; 180 kDa Sin3-associated polypeptide; Sin3-associated polypeptide p180; Breast cancer-associated antigen BRCAA1; Hi
Human Swiss Prot No.:	Q4LE39
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BRCAA1 Polyclonal Antibody detects endogenous levels of BRCAA1 protein.
Storage and Stability:	-20°C/1 year