

Antibodies Product Datasheet

Product Name:	Brg-1 Antibody
Catalog Number:	R12-2044
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Brg-1. AA range:1565-1614
Concentration:	1 mg/ml
Gene Name:	SMARCA4
Protein Name:	Transcription activator BRG1
Other Name:	SMARCA4; BAF190A; BRG1; SNF2B; SNF2L4; Transcription activator BRG1; ATP-dependent helicase SMARCA4; BRG1-associated factor 190A; BAF190A; Mitotic growth and transcription activator; Protein BRG-1; Protein brahma homolog 1; SNF2-beta; SWI/S
Human Swiss Prot No.:	P51532
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Brg-1 Polyclonal Antibody detects endogenous levels of Brg-1 protein.
Storage and Stability:	-20°C/1 year