

Antibodies Product Datasheet

Product Name:	Brm Antibody
Catalog Number:	R12-2045
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Brm. AA range:1328-1377
Concentration:	1 mg/ml
Gene Name:	SMARCA2
Protein Name:	Probable global transcription activator SNF2L2
Other Name:	SMARCA2; BAF190B; BRM; SNF2A; SNF2L2; Probable global transcription activator SNF2L2; ATP-dependent helicase SMARCA2; BRG1-associated factor 190B; BAF190B;Protein brahma homolog; hBRM; SNF2-alpha; SWI/SNF-related matrix-associated actin-dependent regulato
Human Swiss Prot No.:	P51531
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Brm Polyclonal Antibody detects endogenous levels of Brm protein.
Storage and Stability:	-20°C/1 year