

Antibodies Product Datasheet

Product Name:	DFFB Antibody
Catalog Number:	R12-2049
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DFFB
Concentration:	1 mg/ml
Gene Name:	DFFB CAD DFF2 DFF40
Protein Name:	DFFB
Other Name:	DNA fragmentation factor subunit beta (EC 3.-.-) (Caspase-activated deoxyribonuclease) (CAD) (Caspase-activated DNase) (Caspase-activated nuclease) (CPAN) (DNA fragmentation factor 40 kDa subunit) (DFF-40)
Human Swiss Prot No.:	O76075
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DFFB
Storage and Stability:	-20°C/1 year