

Antibodies Product Datasheet

Product Name:	CD179b Antibody
Catalog Number:	R12-2066
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD179b. AA range:26-75
Concentration:	1 mg/ml
Gene Name:	IGLL1/IGLC1/IGLC2/IGLC3/IGLC6/IGLC7
Protein Name:	Immunoglobulin lambda-like polypeptide 1
Other Name:	IGLL1; IGL1; Immunoglobulin lambda-like polypeptide 1; CD179 antigen-like family member B; Ig lambda-5; Immunoglobulin omega polypeptide; Immunoglobulin-related protein 14.1; CD antigen CD179b; IGLC1; Ig lambda-1 chain C regions; IGLC2; Ig
Human Swiss Prot No.:	P15814/P0CG04/P0CG05/P0CG06/P0CF74/A0M8Q6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD179b Polyclonal Antibody detects endogenous levels of CD179b protein.
Storage and Stability:	-20°C/1 year