

Antibodies Product Datasheet

Product Name:	CD34 Antibody
Catalog Number:	R12-2075
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD34. AA range:208-257
Concentration:	1 mg/ml
Gene Name:	CD34
Protein Name:	Hematopoietic progenitor cell antigen CD34
Other Name:	CD34; Hematopoietic progenitor cell antigen CD34; CD antigen CD34
Human Swiss Prot No.:	P28906
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD34 Polyclonal Antibody detects endogenous levels of CD34 protein.
Storage and Stability:	-20°C/1 year