

Antibodies Product Datasheet

Product Name:	CHD1L Antibody
Catalog Number:	R12-2091
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CHD1L. AA range:545-594
Concentration:	1 mg/ml
Gene Name:	CHD1L
Protein Name:	Chromodomain-helicase-DNA-binding protein 1-like
Other Name:	CHD1L; ALC1; Chromodomain-helicase-DNA-binding protein 1-like; Amplified in liver cancer protein 1
Human Swiss Prot No.:	Q86WJ1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CHD1L Polyclonal Antibody detects endogenous levels of CHD1L protein.
Storage and Stability:	-20°C/1 year