

Antibodies Product Datasheet

Product Name:	Connexin-32 Antibody
Catalog Number:	R12-2099
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Connexin-32. AA range:66-115
Concentration:	1 mg/ml
Gene Name:	GJB1
Protein Name:	Gap junction beta-1 protein
Other Name:	GJB1; CX32; Gap junction beta-1 protein; Connexin-32; Cx32; GAP junction 28 kDa liver protein
Human Swiss Prot No.:	P08034
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Connexin-32 Polyclonal Antibody detects endogenous levels of Connexin-32 protein.
Storage and Stability:	-20°C/1 year