

Antibodies Product Datasheet

Product Name:	CRF Antibody
Catalog Number:	R12-2103
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CRF. AA range:104-153
Concentration:	1 mg/ml
Gene Name:	CRH
Protein Name:	Corticoliberin
Other Name:	CRH; Corticoliberin; Corticotropin-releasing factor; CRF; Corticotropin-releasing hormone
Human Swiss Prot No.:	P06850
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRF Polyclonal Antibody detects endogenous levels of CRF protein.
Storage and Stability:	-20°C/1 year