

Antibodies Product Datasheet

Product Name:	Defensin α 3 Antibody
Catalog Number:	R12-2113
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Defensin alpha3. AA range:3-52
Concentration:	1 mg/ml
Gene Name:	DEFA3
Protein Name:	Neutrophil defensin 3
Other Name:	DEFA3; DEF3; Neutrophil defensin 3; Defensin; alpha 3; HNP-3; HP-3; HP3
Human Swiss Prot No.:	P59666
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Defensin α 3 Polyclonal Antibody detects endogenous levels of Defensin α 3 protein.
Storage and Stability:	-20°C/1 year