

Antibodies Product Datasheet

Product Name:	DP-2 Antibody
Catalog Number:	R12-2118
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DP-2. AA range:64-113
Concentration:	1 mg/ml
Gene Name:	TFDP2
Protein Name:	Transcription factor Dp-2
Other Name:	TFDP2; DP2; Transcription factor Dp-2; E2F dimerization partner 2
Human Swiss Prot No.:	Q14188
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DP-2 Polyclonal Antibody detects endogenous levels of DP-2 protein.
Storage and Stability:	-20°C/1 year