

Antibodies Product Datasheet

Product Name:	ECP Antibody
Catalog Number:	R12-2123
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ECP. AA range:102-151
Concentration:	1 mg/ml
Gene Name:	RNASE3
Protein Name:	Eosinophil cationic protein
Other Name:	RNASE3; ECP; RNS3; Eosinophil cationic protein; ECP; Ribonuclease 3; RNase 3
Human Swiss Prot No.:	P12724
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ECP Polyclonal Antibody detects endogenous levels of ECP protein.
Storage and Stability:	-20°C/1 year