

Antibodies Product Datasheet

Product Name:	EEA1 Antibody
Catalog Number:	R12-2124
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EEA1. AA range:205-254
Concentration:	1 mg/ml
Gene Name:	EEA1
Protein Name:	Early endosome antigen 1
Other Name:	EEA1; ZFYVE2; Early endosome antigen 1; Endosome-associated protein p162; Zinc finger FYVE domain-containing protein 2
Human Swiss Prot No.:	Q15075
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EEA1 Polyclonal Antibody detects endogenous levels of EEA1 protein.
Storage and Stability:	-20°C/1 year