

Antibodies Product Datasheet

Product Name:	eIF5A2 Antibody
Catalog Number:	R12-2128
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human eIF5A2. AA range:78-127
Concentration:	1 mg/ml
Gene Name:	EIF5A2
Protein Name:	Eukaryotic translation initiation factor 5A-2
Other Name:	EIF5A2; Eukaryotic translation initiation factor 5A-2; eIF-5A-2; eIF-5A2; Eukaryotic initiation factor 5A isoform 2
Human Swiss Prot No.:	Q9GZV4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF5A2 Polyclonal Antibody detects endogenous levels of eIF5A2 protein.
Storage and Stability:	-20°C/1 year