

Antibodies Product Datasheet

Product Name:	GABAA R β 2 Antibody
Catalog Number:	R12-2152
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GABAA R β 2. AA range:382-431
Concentration:	1 mg/ml
Gene Name:	GABRB2
Protein Name:	Gamma-aminobutyric acid receptor subunit beta-2
Other Name:	GABRB2; Gamma-aminobutyric acid receptor subunit beta-2; GABA(A) receptor subunit beta-2
Human Swiss Prot No.:	P47870
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABAA R β 2 Polyclonal Antibody detects endogenous levels of GABAA R β 2 protein.
Storage and Stability:	-20°C/1 year