

Antibodies Product Datasheet

Product Name:	GnRH I Antibody
Catalog Number:	R12-2167
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000 ,Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GnRH I. AA range:43-92
Concentration:	1 mg/ml
Gene Name:	GNRH1
Protein Name:	Progonadoliberin-1
Other Name:	GNRH1; GNRH; GRH; LHRH; Progonadoliberin-1; Progonadoliberin I
Human Swiss Prot No.:	P01148
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GnRH I Polyclonal Antibody detects endogenous levels of GnRH I protein.
Storage and Stability:	-20°C/1 year