

Antibodies Product Datasheet

Product Name:	GRO γ Antibody
Catalog Number:	R12-2172
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GRO γ . AA range:58-107
Concentration:	1 mg/ml
Gene Name:	CXCL3
Protein Name:	C-X-C motif chemokine 3
Other Name:	CXCL3; GRO3; GROG; SCYB3; C-X-C motif chemokine 3; GRO-gamma(1-73); Growth-regulated protein gamma; GRO-gamma; Macrophage inflammatory protein 2-beta; MIP2-beta
Human Swiss Prot No.:	P19876
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GRO γ Polyclonal Antibody detects endogenous levels of GRO γ protein.
Storage and Stability:	-20°C/1 year