

Antibodies Product Datasheet

Product Name:	GRP1 Antibody
Catalog Number:	R12-2173
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GRP1. AA range:351-400
Concentration:	1 mg/ml
Gene Name:	CYTH3
Protein Name:	Cytohesin-3
Other Name:	CYTH3; ARNO3; GRP1; PSCD3; Cytohesin-3; ARF nucleotide-binding site opener 3; Protein ARNO3; General receptor of phosphoinositides 1; Grp1; PH; SEC7 and coiled-coil domain-containing protein 3
Human Swiss Prot No.:	O43739
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GRP1 Polyclonal Antibody detects endogenous levels of GRP1 protein.
Storage and Stability:	-20°C/1 year