

Antibodies Product Datasheet

Product Name:	HLA Class I Antibody
Catalog Number:	R12-2182
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HLA Class I. AA range:204-253
Concentration:	1 mg/ml
Gene Name:	HLA-A HLAA
Protein Name:	major histocompatibility complex, class I
Other Name:	HLA-A HLAA
Human Swiss Prot No.:	P01891/P01892/P04439/P05534/P10314/P10316/P13746/P16189/P18462/P30443/P30450/P30453/P30455/P30456/P30457/P30512/Q09160
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous HLA Class I protein
Storage and Stability:	-20°C/1 year