

Antibodies Product Datasheet

Product Name:	HoxB5 Antibody
Catalog Number:	R12-2186
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HoxB5. AA range:150-199
Concentration:	1 mg/ml
Gene Name:	HOXB5
Protein Name:	Homeobox protein Hox-B5
Other Name:	HOXB5; HOX2A; Homeobox protein Hox-B5; Homeobox protein HHO.C10; Homeobox protein Hox-2A; Homeobox protein Hu-1
Human Swiss Prot No.:	P09067
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HoxB5 Polyclonal Antibody detects endogenous levels of HoxB5 protein.
Storage and Stability:	-20°C/1 year