

Antibodies Product Datasheet

Product Name:	HSC 70 Antibody
Catalog Number:	R12-2188
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HSC 70. AA range:202-251
Concentration:	1 mg/ml
Gene Name:	HSPA8
Protein Name:	Heat shock cognate 71 kDa protein
Other Name:	HSPA8; HSC70; HSP73; HSPA10; Heat shock cognate 71 kDa protein; Heat shock 70 kDa protein 8
Human Swiss Prot No.:	P11142
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HSC 70 Polyclonal Antibody detects endogenous levels of HSC 70 protein.
Storage and Stability:	-20°C/1 year