

Antibodies Product Datasheet

Product Name:	IDHP Antibody
Catalog Number:	R12-2190
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IDHP
Concentration:	1 mg/ml
Gene Name:	IDH2
Protein Name:	IDHP
Other Name:	Isocitrate dehydrogenase [NADP], mitochondrial (IDH) (EC 1.1.1.42) (ICD-M) (IDP) (NADP(+)-specific ICDH) (Oxalosuccinate decarboxylase)
Human Swiss Prot No.:	P48735
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human IDHP
Storage and Stability:	-20°C/1 year