

Antibodies Product Datasheet

Product Name:	Ihh Antibody
Catalog Number:	R12-2198
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Ihh. AA range:209-258
Concentration:	1 mg/ml
Gene Name:	IHH
Protein Name:	Indian hedgehog protein
Other Name:	IHH; Indian hedgehog protein; IHH; HHG-2
Human Swiss Prot No.:	Q14623
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ihh Polyclonal Antibody detects endogenous levels of Ihh protein.
Storage and Stability:	-20°C/1 year