

Antibodies Product Datasheet

Product Name:	IL-32 Antibody
Catalog Number:	R12-2209
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IL-32. AA range:177-226
Concentration:	1 mg/ml
Gene Name:	IL32
Protein Name:	Interleukin-32
Other Name:	IL32; NK4; TAIF; Interleukin-32; IL-32; Natural killer cells protein 4; Tumor necrosis factor alpha-inducing factor
Human Swiss Prot No.:	P24001
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-32 Polyclonal Antibody detects endogenous levels of IL-32 protein.
Storage and Stability:	-20°C/1 year