

Antibodies Product Datasheet

Product Name:	ING4 Antibody
Catalog Number:	R12-2215
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ING4. AA range:107-156
Concentration:	1 mg/ml
Gene Name:	ING4
Protein Name:	Inhibitor of growth protein 4
Other Name:	ING4; My036; Inhibitor of growth protein 4; p29ING4
Human Swiss Prot No.:	Q9UNL4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ING4 Polyclonal Antibody detects endogenous levels of ING4 protein.
Storage and Stability:	-20°C/1 year