

Antibodies Product Datasheet

Product Name:	Integrin α 3 Antibody
Catalog Number:	R12-2220
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Integrin alpha3. AA range:482-531
Concentration:	1 mg/ml
Gene Name:	ITGA3
Protein Name:	Integrin alpha-3
Other Name:	ITGA3; MSK18; Integrin alpha-3; CD49 antigen-like family member C; FRP-2; Galactoprotein B3; GAPB3; VLA-3 subunit alpha; CD antigen CD49c
Human Swiss Prot No.:	P26006
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Integrin α 3 Polyclonal Antibody detects endogenous levels of Integrin α 3 protein.
Storage and Stability:	-20°C/1 year