

Antibodies Product Datasheet

Product Name:	MIP-1 α Antibody
Catalog Number:	R12-2239
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MIP-1 α . AA range:26-75
Concentration:	1 mg/ml
Gene Name:	CCL3
Protein Name:	C-C motif chemokine 3
Other Name:	CCL3; G0S19-1; MIP1A; SCYA3; C-C motif chemokine 3; G0/G1 switch regulatory protein 19-1; Macrophage inflammatory protein 1- α ; MIP-1- α ; PAT 464.1; SIS-beta; Small-inducible cytokine A3; Tonsillar lymphocyte LD78 α protein
Human Swiss Prot No.:	P10147
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MIP-1 α Polyclonal Antibody detects endogenous levels of MIP-1 α protein.
Storage and Stability:	-20°C/1 year