

Antibodies Product Datasheet

Product Name:	ML-IAP Antibody
Catalog Number:	R12-2241
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ML-IAP. AA range:162-211
Concentration:	1 mg/ml
Gene Name:	BIRC7
Protein Name:	Baculoviral IAP repeat-containing protein 7
Other Name:	BIRC7; KIAP; LIVIN; MLIAP; RNF50; Baculoviral IAP repeat-containing protein 7; Kidney inhibitor of apoptosis protein; KIAP; Livin; Melanoma inhibitor of apoptosis protein; ML-IAP; RING finger protein 50
Human Swiss Prot No.:	Q96CA5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ML-IAP Polyclonal Antibody detects endogenous levels of ML-IAP protein.
Storage and Stability:	-20°C/1 year