

Antibodies Product Datasheet

Product Name:	NET-5 Antibody
Catalog Number:	R12-2255
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NET-5. AA range:72-121
Concentration:	1 mg/ml
Gene Name:	TSPAN9
Protein Name:	Tetraspanin-9
Other Name:	TSPAN9; NET5; Tetraspanin-9; Tspan-9; Tetraspan NET-5
Human Swiss Prot No.:	O75954
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NET-5 Polyclonal Antibody detects endogenous levels of NET-5 protein.
Storage and Stability:	-20°C/1 year