

Antibodies Product Datasheet

Product Name:	NF90 Antibody
Catalog Number:	R12-2260
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NF90. AA range:302-351
Concentration:	1 mg/ml
Gene Name:	ILF3
Protein Name:	Interleukin enhancer-binding factor 3
Other Name:	ILF3; DRBF; MPHOSPH4; NF90; Interleukin enhancer-binding factor 3; Double-stranded RNA-binding protein 76; DRBP76; M-phase phosphoprotein 4; MPP4;Nuclear factor associated with dsRNA; NFAR; Nuclear factor of activated T-cells 90 kDa; NF-AT-90; Translation
Human Swiss Prot No.:	Q12906
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NF90 Polyclonal Antibody detects endogenous levels of NF90 protein.
Storage and Stability:	-20°C/1 year