

Antibodies Product Datasheet

Product Name:	RAP1 Antibody
Catalog Number:	R12-2330
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RAP1. AA range:203-252
Concentration:	1 mg/ml
Gene Name:	TERF2IP
Protein Name:	Telomeric repeat-binding factor 2-interacting protein 1
Other Name:	TERF2IP; DRIP5; RAP1; PP8000; Telomeric repeat-binding factor 2-interacting protein 1; TERF2-interacting telomeric protein 1; TRF2-interacting telomeric protein 1; Dopamine receptor-interacting protein 5; Repressor/activator protein 1 homol
Human Swiss Prot No.:	Q9NYB0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RAP1 Polyclonal Antibody detects endogenous levels of RAP1 protein.
Storage and Stability:	-20°C/1 year