

Antibodies Product Datasheet

Product Name:	SDHA Antibody
Catalog Number:	R12-2348
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SDHA. AA range:551-600
Concentration:	1 mg/ml
Gene Name:	SDHA
Protein Name:	Succinate dehydrogenase [ubiquinone] flavoprotein subunit mitochondrial
Other Name:	SDHA; SDH2; SDHF; Succinate dehydrogenase [ubiquinone] flavoprotein subunit; mitochondrial; Flavoprotein subunit of complex II; Fp
Human Swiss Prot No.:	P31040
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SDHA Polyclonal Antibody detects endogenous levels of SDHA protein.
Storage and Stability:	-20°C/1 year