

Antibodies Product Datasheet

Product Name:	SP-B Antibody
Catalog Number:	R12-2366
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SP-B. AA range:243-292
Concentration:	1 mg/ml
Gene Name:	SFTPB
Protein Name:	Pulmonary surfactant-associated protein B
Other Name:	SFTPB; SFTP3; Pulmonary surfactant-associated protein B; SP-B; 18 kDa pulmonary-surfactant protein; 6 kDa protein; Pulmonary surfactant-associated proteolipid SPL(Phe)
Human Swiss Prot No.:	P07988
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SP-B Polyclonal Antibody detects endogenous levels of SP-B protein.
Storage and Stability:	-20°C/1 year