

Antibodies Product Datasheet

Product Name:	VAMP-1/2/3 Antibody
Catalog Number:	R12-2405
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human VAMP-1/2/3. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	VAMP1/VAMP2/VAMP3
Protein Name:	Vesicle-associated membrane protein 1/2/3
Other Name:	VAMP1; SYB1; Vesicle-associated membrane protein 1; VAMP-1; Synaptobrevin-1; VAMP2; SYB2; Vesicle-associated membrane protein 2; VAMP-2; Synaptobrevin-2; VAMP3; SYB3; Vesicle-associated membrane protein 3; VAMP-3; Cellubrevin; CEB; Synaptob
Human Swiss Prot No.:	P23763/P63027/Q15836
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	VAMP-1/2/3 Polyclonal Antibody detects endogenous levels of VAMP-1/2/3 protein.
Storage and Stability:	-20°C/1 year