
Antibodies Product Datasheet

Product Name:	ZP4 Antibody
Catalog Number:	R12-3706
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	ZP4 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZP4.
Concentration:	1 mg/ml
Gene Name:	ZP4
Other Name:	Zona pellucida sperm-binding protein 4; Zona pellucida sperm-binding protein B; Contains Processed zona pellucida sperm-binding protein 4; ZP4; ZPB
Human Swiss Prot No.:	Q12836
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ZP4 Antibodyantibody detects endogenous levels of ZP4.
Storage and Stability:	Stable at -20°C for at least 1 year.