

Antibodies Product Datasheet

Product Name:	FITC* Goat Anti Rabbit IgG(H+L)
Catalog Number:	SA0004
Reactivity:	Rabbit
Conjugate:	FITC
Applications:	WB,IHC,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:500-1000 IF 1:500-2000
Modification:	Not Modification Specific
Source:	Goat
Purification:	The antibody was affinity-purified from goat antiserum by affinity-chromatography using epitope-specific immunogen.
Concentration:	4 mg/ml
Other Name:	FITC* Goat Anti-Rabbit IgG(H+L)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year