

Antibodies Product Datasheet

Product Name:	HRP*Polyclonal Rabbit Anti Chicken IgY (IgG)
Catalog Number:	SA0104
Reactivity:	ChickenIgY(IgG)
Conjugate:	HRP
Applications:	WB,IHC-p,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
Modification:	Not Modification Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Concentration:	4-5mg/mL
Other Name:	HRP*Polyclonal Rabbit Anti-Chicken IgY (IgG)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year