

Antibodies Product Datasheet

Product Name:	HRP*Polyclonal Rabbit Anti Horse IgG(h+l)
Catalog Number:	SA0105
Reactivity:	HorseIgG(h+l)
Conjugate:	HRP
Applications:	WB,IHC-p,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
Modification:	Not Modification Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Concentration:	4-5mg/mL
Other Name:	HRP*Polyclonal Rabbit Anti-Horse IgG(h+l)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year

SA0105 | HRP*Polyclonal Rabbit Anti Horse IgG(h+l)

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022