

Antibodies Product Datasheet

Product Name:	AU*Polyclonal Goat Anti Chicken IgY (IgG)(h+l)
Catalog Number:	SA0182
Reactivity:	ChickenIgY(IgG)(h+l)
Conjugate:	AU
Applications:	WB,IHC-p,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
Modification:	Not Modification Specific
Source:	Goat
Purification:	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Concentration:	4-5mg/mL
Other Name:	AU*Polyclonal Goat Anti-Chicken IgY (IgG)(h+l)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year