

Antibodies Product Datasheet

Product Name:	FITC*Polyclonal Goat Anti Human IgM(u chain)
Catalog Number:	SA0245
Reactivity:	HumanIgM(uchain)
Conjugate:	FITC
Applications:	WB,IHC-p,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
Modification:	Not Modification Specific
Source:	Goat
Purification:	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Concentration:	4-5mg/mL
Other Name:	FITC*Polyclonal Goat Anti-Human IgM(u chain)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year