

Antibodies Product Datasheet

Product Name:	BIOTIN*Polyclonal Rabbit Anti Donkey IgG (h+I)
Catalog Number:	SA0307
Reactivity:	DonkeyIgG(h+I)
Conjugate:	Biotin
Applications:	WB,IHC-p,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
Modification:	Not Modification Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Concentration:	4-5mg/mL
Other Name:	BIOTIN*Polyclonal Rabbit Anti- Donkey IgG (h+I)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year